IgM-Reducing Assay Diluent

Minimizes sample matrix differentials and IgM-associated false positive signal.

IgM-Reducing Assay Diluent is formulated to ensure a more accurate and consistent ELISA performance when testing serum and plasma samples in all sandwich ELISA formats. This novel assay diluent utilizes special reducing agents to minimize IgM-associated false positive signal. Rogue IgM in samples can mediate non-specific bridging of the HRP conjugate to coated and

Equalize the matrix with an assay diluent diluent to contains the antibody small molecules every well to target antigen equalize the and interfering binds to in the sample the target may suppress molecules that antibody to antigen antigen binding mix with the assay diluent efficiency

blocked ELISA plate surfaces, creating confusing false positive signal that can easily invalidate any ELISA-derived sample data.

IgM-Reducing Assay Diluent provides a mechanism for equalization of matrix complexity differences between the standard curve diluent and test samples, which would otherwise result in under-recovery of the target analyte present in sample wells. Complex and concentrated protein environments (the matrix) present in serum or plasma samples, unlike the simple matrix of the standards, will greatly reduce the antigen binding efficiency of the plate-adsorbed antibodies, resulting in a gross underestimation of the amount of target analyte present in the test samples. This serious issue common to serum and plasma ELISA format assays may be addressed with the use of an appropriate assay diluent. IgM-Reducing Assay Diluent also contains additives to inhibit complement and thrombin activity present in serum and plasma samples.

To use, simply add 50-100 μ L to every well of the ELISA plate, including all wells designated for standards, controls, and samples. Then add the standards, controls, and samples to the plate. An antimicrobial agent is included in the product formulation to allow for room temperature bench top use and extensive storage stability at 2-8°C.

IGM-REDUCING ASSAY DILUENT

Catalog#
LS-M20-100
LS-M20-500
LS-M20-1

INSTRUCTIONS:

- Dilute the standard curve, controls, and the samples as necessary. LSBio offers several formulations of sample diluents in which to prepare the samples.
- 2. Pipette 50-100 μL IgM-Reducing Assay Diluent per well into every well of the plate.
- 3. Pipette 50-200 µL of each standard, control, and sample into the plate.
- 4. Run the rest of the assay according to the specific ELISA protocol.
- 5. Analyze the data. Because all the wells, including the standards and controls, received the same volume of assay diluent, there is no need to account for this dilution when calculating the results.

For more ELISA protocols and information, please visit www.LSBio.com

SPECIFICATIONS:

- · Clear to cloudy, slightly white liquid
- 1X ready to use
- pH 6.4-6.6

STORAGE:

- Shelf life is 18 months at 2-8°C
- 1 week at room temperature

SAFETY & USAGE:

- Danger!
- Contains ≥ 3 < 5% sodium metabisulfite
- Causes serious eye damage
- Product also includes ≤ 0.1% sodium azide
- SDS available upon request
- Not for human or drug use
- For research use only



2401 4th Avenue, Suite 900, Seattle, WA 98121 **USA TOLL-FREE** 1-866-819-4732 **LOCAL** 206-374-1102