

# Nitric Oxide (NO) Synthase Assay Kit (Colorimetric)

LS-K196-100 (100 Tests) • Store at -20°C



## Introduction

Nitric oxide (NO) is a reactive radical that plays an important role in many key physiological functions. NO, an oxidation product of arginine by nitric oxide synthase (NOS), is involved in host defense and development, activation of regulatory proteins and direct covalent interaction with functional biomolecules. Simple, direct and non-radioactive procedures for measuring NOS are becoming popular in Research and Drug Discovery. This Nitric Oxide Synthase Assay Kit involves two steps: a NOS reaction step during which NO is produced followed by an NO detection step. Since the NO generated by NOS is rapidly oxidized to nitrite and nitrate, the NO production is measured following reduction of nitrate to nitrite using an improved Griess Method. The procedure is simple and the time required for sample pretreatment and assay is reduced to as short as 40 min.

## Key Features

- Sensitive and accurate. Detection range 0.25 - 25 U/L in 96-well plate.
- Rapid and reliable. Can be completed in 40 min if reduction of  $\text{NO}_3^-$  to  $\text{NO}_2^-$  is performed at 60°C.

## Applications

- Direct Assays: NOS activity in biological samples.
- Drug Discovery/Pharmacology: effects of drugs on NOS activity.

## Components

Component	K196-100
	100 Tests
Assay Buffer	10 mL
Substrate	600 $\mu\text{L}$
GDH	120 $\mu\text{L}$
Reagent A	12 mL
Reagent B	500 $\mu\text{L}$
Reagent C	12 mL
Reagent D (2 Tubes)	Dried
Reagent E	1.5 mL
$\text{ZnSO}_4$	1 mL
Standard	1 mL
NaOH	1 mL

## Materials Not Supplied

Pipetting devices, Eppendorf tubes, Eppendorf centrifuge, clear, flat-bottomed 96 well plates or cuvettes, plate reader or spectrophotometer and heat block or hot water bath (optional).

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)

# Nitric Oxide (NO) Synthase Assay Kit (Colorimetric)

LS-K196-100 (100 Tests) • Store at -20°C



## Storage

The kit is shipped on ice. Store Assay Buffer, Substrate, Reagent D, Reagent E and GDH at -20°C. Store all other reagents at 4°C. Shelf life of six months after receipt. Use Reagent D within 1 week after reconstitution.

## Assay Procedure

Prior to assay, equilibrate all components to room temperature. Reconstitute one tube of Reagent D with 300 µL dH<sub>2</sub>O. (If assaying more than 60 wells, reconstitute both tubes of Reagent D) Store unused reconstituted Reagent D at -20°C and use within 1 week. Prewarm Assay Buffer to 37°C. Keep GDH on ice. If precipitates are present in Reagent B, warm at 37°C until redissolved (~10-15 min).

## Sample treatment

Tissue or cell samples are homogenized in 1 x PBS (pH 7.4). Centrifuge at 10,000g or higher at 4°C. Use supernatant for NOS assay.

## Standard preparation

Prepare 200 µL 500 µM Premix by mixing 100 µL 1.0 mM Standard and 100 µL distilled water. Dilute standards in 1.5 mL centrifuge tubes as described in the Table.

No	Premix + H <sub>2</sub> O	Nitrite (µM)
1	50 µL + 0 µL	500
2	30 µL + 20 µL	300
3	15 µL + 35 µL	150
4	0 µL + 50 µL	0

## Procedure

**NOS Reaction:** If samples will not require deproteinization (i.e. purified NOS), add 20 µL of each sample and standard to separate labeled Eppendorf tubes. Each sample requires at least two tubes: one reaction tube and one sample blank tube. Immediately prior to starting the reaction, prepare enough NOS Working Reagent (NOS WR) for all sample reaction tubes and standards by mixing per reaction tube: 65 µL Assay Buffer, 4 µL Substrate, 4 µL Reconstituted Reagent D, 10 µL Reagent E and 1 µL GDH. For the sample blanks, use 8 µL dH<sub>2</sub>O instead of the Substrate and Reagent D. Add 80 µL of the appropriate NOS WR to each tube and incubate at 37°C for 20 min. After 20 min immediately add 200 µL of the NO Detection Reagent (NO DR) (see next section: NO Measurement) to each tube to kill the NOS reaction.

For samples requiring deproteinization which include serum, plasma, whole blood, cell culture media containing FBS, tissue or cell lysates, add 25 µL of each sample and standard to separate labeled Eppendorf tubes. Each sample requires at least two tubes: one reaction tube and one sample blank tube. Immediately prior to starting the reaction, prepare enough NOS WR for all sample reaction tubes and standards by mixing per reaction tube: 80 µL Assay Buffer, 5 µL Substrate, 5 µL Reconstituted Reagent D, 13 µL Reagent E and 1 µL GDH. For the sample blanks, use 10 µL dH<sub>2</sub>O instead of the Substrate and Reagent D. Add 100 µL of the appropriate NOS WR to each tube and incubate at 37°C for 20 min. After 20 min immediately proceed to the deproteinization step.

**Deproteinization:** Add 7 µL ZnSO<sub>4</sub> to each sample and standard tube. Vortex and then add 7 µL NaOH. Vortex again and centrifuge 10 min at 14,000 rpm. Transfer 100 µL of the clear supernatant to a clean tube and proceed to the NO Measurement step.

**NO Measurement:** Immediately prior to starting the reaction, prepare enough NO Detection Reagent (NO DR) for all samples and standards by mixing per reaction tube: 100 µL Reagent A, 4 µL Reagent B and 100 µL Reagent C. Add 200 µL of the NO DR to each sample and standard tube and incubate for 5 min at 60°C. (Alternatively, the reaction can be run at 37°C for 60 min or RT for 150 min.)

Briefly centrifuge the reaction tubes to pellet any condensation and transfer 250 µL of each reaction to separate wells in a 96 well plate. Read OD at 500-570nm (peak 540 nm).

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)

# Nitric Oxide (NO) Synthase Assay Kit (Colorimetric)

LS-K196-100 (100 Tests) • Store at -20°C



---

## Calculations

Subtract blank OD (Std 4) from the standard OD values and plot the OD against standard concentrations. Determine the slope using linear regression fitting. The NOS activity of the Sample is then calculated as

$$\text{NOS Activity} = \frac{\text{OD}_{\text{SAMPLE}} - \text{OD}_{\text{BLANK}}}{\text{Slope}} \times \frac{1}{t} \text{ (U/L)}$$

OD<sub>SAMPLE</sub> and OD<sub>BLANK</sub> are optical density values of the sample and sample blank, respectively. *t* is the reaction time (20 min).

Unit definition: one unit of NOS catalyzes the production of 1 μmole of nitric oxide per minute under the assay conditions (pH 7.5 and 37°C).

## General Considerations

Antioxidants and nucleophiles (e.g. β-mercaptoethanol, glutathione, dithiothreitol and cysteine) may interfere with this assay. Avoid using these compounds during sample preparation. However, if β-mercaptoethanol or dithiothreitol must be used, an equal concentration needs to be added to the standards.

Version: V.08.09.2018

---

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)