# Catalase Assay Kit (Colorimetric/Fluorometric)

LS-K245-100 (100 Tests) • Store at -20°C



### Introduction

CATALASE (EC 1.11.1.6), is a ubiquitous antioxidant enzyme that catalyzes the decomposition of hydrogen peroxide  $(H_2O_2)$  to water and oxygen.



By preventing excessive  $H_2O_2$  build up, catalase allows important cellular processes which produce  $H_2O_2$  as a byproduct to occur safely. Deficiency in catalase activity has been associated with grey hair and peroxisomal disorder acatalasia. Simple, direct and high-throughput assays for catalase activity find wide applications. LSBio's' improved assay directly measures catalase degradation of  $H_2O_2$  using a redox dye. The change in color intensity at 570nm or fluorescence intensity ( $\lambda$ em/ex = 585/530nm) is directly proportional to the catalase activity in the sample.

### **Key Features**

- Sensitive and accurate. Use 10 µL sample. Linear detection range 0.2 to 5 U/L catalase activity.
- Simple and Convenient. The procedure involves adding a Substrate to sample, incubation for 30 min, followed by a Detection Reagent and reading the optical density or fluorescence intensity.

### Applications

- Direct Assays: catalase activity in biological samples e.g. serum, plasma, urine, saliva, cell culture etc.
- Drug Discovery/Pharmacology: effects of drugs on catalase activity.

	K245-100	
Component	100 Tests	
Assay Buffer	25 mL	
HRP Enzyme	120 μL	
Dye Reagent	120 μL	
$H_2O_2$ Solution	$100 \ \mu L \ 3\% \ H_2O_2$	
Positive Control	8 μL Catalase	

#### Components

# **Materials Not Supplied**

Pipetting devices, flat-bottom 96-well plates, plate reader.

#### Storage

The kit is shipped on ice. Store all kit components at -20 °C.

# Catalase Assay Kit (Colorimetric/Fluorometric)

LS-K245-100 (100 Tests) • Store at -20°C



# **Assay Procedure**

1. *Reagent Preparation*. Equilibrate all components to room temperature. Briefly centrifuge all tubes before opening. Keep thawed HRP Enzyme on ice.

For colorimetric assays, use a clear flat-bottom 96-well plate. For fluorometric assays, use a solid black flat-bottom 96-well plate.

Samples and Controls: transfer 10 µL sample into wells of the 96-well plate. In addition, for each assay run, prepare one sample blank well that contains only 10 µL Assay Buffer.

Add 400  $\mu$ L Assay Buffer to Positive Control tube and mix well. Transfer 10  $\mu$ L of the reconstituted Positive Control into separate wells.

Note: (1). For unknown samples, perform several dilutions to ensure that catalase activity is within the linear range 0.2 to 5 U/L. (2) The provided catalase serves as a positive control to ensure assay is working and should not be used to calculate the Sample catalase activity.

- Enzyme Reaction. Mix 5 μL 3% H2O2 and 914 μL dH<sub>2</sub>O (final 4.8 mM). Prepare enough 50 μM H<sub>2</sub>O<sub>2</sub> Substrate for sample, positive control and sample blank by mixing, for each well, 1 μL of the 4.8 mM H<sub>2</sub>O<sub>2</sub> with95 μL Assay Buffer. Note: diluted H<sub>2</sub>O<sub>2</sub> is not stable. Prepare fresh dilutions for each experiment. Add 90 μL of the 50 μM Substrate to these wells to initiate the catalase reaction. Tap plate quick to mix. Incubate 30 min at room temperature. During the incubation time, proceed with Steps 3 and 4 below.
- H<sub>2</sub>O<sub>2</sub> Standard Curve. Mix 40μL of the 4.8 mM H<sub>2</sub>O<sub>2</sub> with 440 μL dH<sub>2</sub>O to yield 400 μM H<sub>2</sub>O<sub>2</sub>. Prepare standards as shown in the Table below. Transfer 10 μL standards into separate wells of the 96-well plate. Add 90 μL Assay Buffer to the standards.

No	400 μM H <sub>2</sub> O <sub>2</sub> + H <sub>2</sub> O	Vol (μL)	H <sub>2</sub> O <sub>2</sub> (μM)
1	100µL +0µL	100	400
2	60µL + 40µL	100	240
3	30µL + 70µL	100	120
4	0μL + 100μL	100	0

- Detection. Prepare enough Detection Reagent by mixing, for each reaction well (Sample, Control and Standard wells), 102 μL Assay Buffer, 1 μL Dye Reagent and 1 μL HRP Enzyme. At the end of the 30 min incubation (Step 2), add 100 μL Detection Reagent per well. Tap plate to mix. Incubate for 10 min.
- 5. Read optical density at 570nm (550 to 585nm) or fluorescence intensity at  $\lambda$ em/ex = 585/530nm.

#### CalculationS

Subtract blank value (#4) from the standard values and plot the  $\Delta$ OD or  $\Delta$ F against standard concentrations. Determine the slope and calculate the catalase activity of Sample,

Catalase (U/L) = 
$$\frac{R_{\text{Sample Blank}} - R_{\text{Sample}}}{\text{Slope } (\mu M^{-1}) \times 30 \text{ min}} \times n$$

R<sub>SAMPLE Blank</sub> and R<sub>SAMPLE</sub> are optical density or fluorescence intensity readings of the Sample Blank and Sample, respectively. Slope is determined from the standard curve. 30 min is the catalase reaction time. *n* is the sample dilution factor.

Unit definition: one unit is the amount of catalase that decomposes  $1 \mu$ mole of H<sub>2</sub>O<sub>2</sub> per min at pH 7.0 and room temperature.

# Catalase Assay Kit (Colorimetric/Fluorometric)



LS-K245-100 (100 Tests) • Store at -20°C

# Sample Data



Version: V.08.09.2018

**FOR RESEARCH USE ONLY!** Not for use in humans. LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121 <u>www.LSBio.com</u> • (206) 464-1554 • <u>TechnicalSupport@LSBio.com</u>