Glutamate Assay Kit (Colorimetric)

LS-K264-100 (100 Tests) • Store at -20°C



Introduction

Glutamate is an important chemical in general metabolism. It is also a crucial mammalian neurotransmitter that is believed to be involved in a number of neurological and psychiatric disorders such as lateral sclerosis, lathyrism, autism and Alzheimer's disease. Glutamate is also widely used as a flavor enhancer in the food industry.

Simple, direct and automation-ready procedures for measuring glutamate concentration are very desirable. LSBio's glutamate assay kit is based on glutamate dehydrogenase catalyzed oxidation of glutamate, in which the formed NADH reduces a formazan (MTT) Reagent. The intensity of the product color, measured at 565 nm, is proportionate to the glutamate concentration in the sample.

Key Features

- Sensitive and accurate. Detection limit of 50 μM, linearity up to 2.5 mM glutamate in 96-well plate assay.
- Convenient. The procedure involves adding a single working reagent, and reading the optical density at time zero and at 30 min at room temperature. No 37°C heater is needed.
- High-throughput. Can be readily automated as a high-throughput 96- well plate assay for thousands of samples per day.

Applications

- Direct Assays: glutamate in serum, plasma, tissue extracts and food extract samples.
- Drug Discovery/Pharmacology: effects of drugs on glutamate levels.

Components

	K264-100
Component	100 Tests
Assay Buffer	10 mL
Enzyme A	120 μL
Enzyme B	120 μL
NAD Solution	1 mL
MTT Solution	1.5 mL
Standard (100 mM Glutamate)	1 mL

Materials Not Supplied

Pipetting (multi-channel) devices. Clear-bottom 96-well plates and plate reader.

Storage

The kit is shipped on ice. Store all kit components at -20 °C.

Glutamate Assay Kit (Colorimetric)

LS-K264-100 (100 Tests) • Store at -20°C



Assay Procedure

1. Calibration Curve. Prepare 600 μ L 2.5 mM Glutamate Premix by mixing 15 μ L 100 mM Standard and 585 μ L distilled water. Dilute standard as follows. Transfer 20 μ L standards into wells of a clear bottom 96-well plate.

No	Premix + H ₂ O	Vol (μL)	Glutamate (mM)
1	100 μL + 0 μL	100	2.5
2	80 μL + 20 μL	100	2.0
3	60 μL + 40 μL	100	1.5
4	40 μL + 60 μL	100	1.0
5	30 μL + 70 μL	100	0.75
6	20 μL + 80 μL	100	0.5
7	10 μL + 90 μL	100	0.25
8	0 μL + 100 μL	100	0

Samples: add 20 μ L sample per well in separate wells. IMPORTANT: Serum and tissue extract samples require a sample blank.

- 2. Reagent Preparation. Spin the Enzyme tubes briefly before pipetting. For each well of reaction, prepare Working Reagent by mixing 60 μL Assay Buffer, 1 μL Enzyme A, 1 μL Enzyme B, 5 μL NAD and 14 μL MTT. Fresh reconstitution is recommended. Where a sample blank in required, prepare a Blank Working Reagent by mixing 60 μL Assay Buffer, 1 μL Enzyme B, 5 μL NAD and 14 μL MTT (i.e. No Enzyme A). This assay is based on a kinetic reaction. To ensure identical incubation time, addition of Working Reagent to samples should be quick and mixing should be brief but thorough. Use of a multi-channel pipettor is recommended. Assays can be executed at room temperature or 30°C.
- 3. Reaction. Add 80 μ L Working Reagent (or Blank Working Reagent where appropriate) per reaction well quickly. Tap plate to mix briefly and thoroughly.
- 4. Read optical density (OD_0) for time "zero" at 565 nm (520-600 nm) and OD_{30} after a 30-min incubation at room temperature.
- 5. Calculation. Subtract OD_0 from OD_{30} for the standard and sample wells. Next, subtract the ΔOD_{water} (Std 8) from each $\Delta OD_{standard}$ and ΔOD_{sample} to obtain the $\Delta \Delta OD_{sample}$. Where a sample blank was required, subtract the ΔOD_{blank} from ΔOD_{sample} to obtain the $\Delta \Delta OD_{sample}$.) Plot the $\Delta \Delta OD_{standard}$'s and use this standard curve to convert the $\Delta \Delta OD_{sample}$ values to sample glutamate concentration.

[Glutamate] =
$$\frac{\Delta \Delta OD_{SAMPLE}}{Slope}$$
 (mM)

Note: If the sample $\Delta\Delta$ OD values are higher than the $\Delta\Delta$ OD value for the 2.5 mM glutamate standard, dilute sample in distilled water and repeat this assay. Multiply the results by the dilution factor.

Conversions: 1 mM glutamate = 14.6 mg/dL.

General Considerations


- 1. This assay is based on an enzyme-catalyzed kinetic reaction. Addition of Working Reagent should be quick and mixing should be brief but thorough. Use of multi-channel pipettor is recommended.
- 2. The following substances interfere and should be avoided in sample preparation: EDTA (>0.5 mM), ascorbic acid, SDS (>0.2%), sodium azide, NP-40 (>1%) and Tween-20 (>1%).

Glutamate Assay Kit (Colorimetric)

LS-K264-100 (100 Tests) • Store at -20°C



Sample Data



Standard Curve in 96-well plate assay

Version: V.08.09.2018