

# Maltose Assay Kit (Colorimetric/Fluorometric)

LS-K295-100 (100 Tests) • Store at -20°C



---

## Introduction

MALTOSE (C<sub>12</sub>H<sub>22</sub>O<sub>11</sub>) is a disaccharide, composed of two glucose units linked by an alpha bond. It is produced from the hydrolysis of glycogen or starch, serving as a source of energy for plants and animals. Maltose can be found in foods such as grains, and other processed products. This maltose assay provides a simple, one step assay for measuring maltose. In this assay, maltose is converted to two glucoses, which are then oxidized to form a colored product. The color intensity of the product at 570 nm or fluorescence at  $\lambda_{ex/em} = 530/585$  nm is directly proportional to maltose concentration in the sample.

## Key Features

- Fast and sensitive. Use of 10  $\mu$ L sample. Linear detection range 2 to 500  $\mu$ M maltose in for colorimetric assays and 1 to 50  $\mu$ M for fluorometric assays.
- Convenient. The procedure involves adding a single working reagent, and reading the absorbance after 60 minutes. Room temperature assay. No 37°C heater is needed.
- High-throughput. "Add-mix-read" type assay. Can be readily automated as a high-throughput 96-well plate assay for thousands of samples per day.

## Applications

- Direct Assays: Maltose in various biological samples such as serum, urine, food and beverages, etc.

## Components

Component	K295-100
	100 Tests
Assay Buffer	10 mL
Enzyme A	120 $\mu$ L
Enzyme Mix	120 $\mu$ L
Dye Reagent	120 $\mu$ L
Standard (5 mM Maltose)	1 mL

## Materials Not Supplied

Pipetting devices, centrifuge tubes, clear flat-bottom 96-well plates (e.g. VWR cat# 82050-760), and plate reader.

## Storage

The kit is shipped on ice. Store all components at -20°C upon receiving. Shelf life: 6 months after receipt.

---

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121  
[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)

# Maltose Assay Kit (Colorimetric/Fluorometric)

LS-K295-100 (100 Tests) • Store at -20°C



---

## Assay Procedure

### Sample Preparation

Clear and slightly colored samples can be assayed directly. It is prudent to test several dilutions to determine an optimal dilution factor  $n$ .

Serum can be assayed directly, after centrifuging first to remove any particulates. Appropriate dilution in distilled water may be required.

Urine samples can be assayed directly, after centrifuging first to remove any particulates, and requires the use of an internal standard.

### Colorimetric Procedure

1. Standards. Prepare 500  $\mu\text{L}$  of 500  $\mu\text{M}$  Premix by mixing 50  $\mu\text{L}$  of the Standard (5 mM) and 450  $\mu\text{L}$  distilled water. Dilute standards in 1.5-mL centrifuge tubes as described in the Table.

No	Premix + H <sub>2</sub> O	Maltose ( $\mu\text{M}$ )
1	100 $\mu\text{L}$ + 0 $\mu\text{L}$	500
2	60 $\mu\text{L}$ + 40 $\mu\text{L}$	300
3	30 $\mu\text{L}$ + 70 $\mu\text{L}$	150
4	0 $\mu\text{L}$ + 100 $\mu\text{L}$	0

Transfer 10  $\mu\text{L}$  standards into separate wells of a clear, flat-bottom 96-well plate. Transfer 10  $\mu\text{L}$  of each sample into separate wells.

2. Samples. Add 10  $\mu\text{L}$  of each sample to two separate wells in a 96 well plate (each sample requires a sample blank).

If using an internal standard, samples will need three separate reactions: 1) sample plus standard, 2) sample alone and 3) sample blank. For the sample plus standard well, add 5  $\mu\text{L}$  500  $\mu\text{M}$  maltose and 10  $\mu\text{L}$  sample. For the sample and sample blank wells, add 5  $\mu\text{L}$  dH<sub>2</sub>O and 10  $\mu\text{L}$  sample.

3. Prepare sufficient Working Reagent (WR) for all sample and standard wells by mixing, for each well: 95  $\mu\text{L}$  Assay Buffer, 1  $\mu\text{L}$  Enzyme A, 1  $\mu\text{L}$  Enzyme Mix, and 1  $\mu\text{L}$  Dye Reagent. For the sample blanks, prepare a Blank WR (BWR) without Enzyme A.

Add 90  $\mu\text{L}$  WR to the four Standards and the Sample Wells and 90  $\mu\text{L}$  BWR to the Sample Blank Wells. Tap plate to mix briefly and thoroughly. Incubate 60 minutes at room temperature.

4. Read optical density at 570 nm (525-605 nm).

### Fluorometric Procedure

1. Standards. Dilute the standards prepared in Colorimetric Procedure 1:10 in dH<sub>2</sub>O. If an internal standard is used, use 5  $\mu\text{L}$  of 50  $\mu\text{M}$  maltose.
2. Transfer 10  $\mu\text{L}$  standards into separate wells of a black 96-well plate. Transfer 10  $\mu\text{L}$  of each sample into separate wells (2 wells per sample if a standard curve is used, 3 wells per sample if an internal standard is used).
3. Add 90  $\mu\text{L}$  of appropriate WR (see Colorimetric Procedure) to each well. Tap plate to mix briefly and thoroughly. Incubate 60 minutes at room temperature.
4. Read fluorescence at  $\lambda_{\text{ex/em}} = 530/585$  nm.

---

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)

# Maltose Assay Kit (Colorimetric/Fluorometric)

LS-K295-100 (100 Tests) • Store at -20°C

## Calculations

Subtract the blank value (#4) from the standard values and plot the  $\Delta OD$  against standard concentrations. Determine the slope and calculate the maltose concentration of the Sample,

$$[\text{Maltose}] = \frac{R_{\text{SAMPLE}} - R_{\text{BLANK}}}{\text{Slope } (\mu\text{M}^{-1})} \times n \quad (\mu\text{M})$$


If an internal standard was used, the sample maltose concentration is computed as follows:

$$[\text{Maltose}] = \frac{R_{\text{SAMPLE}} - R_{\text{BLANK}}}{R_{\text{STANDARD}} - R_{\text{SAMPLE}}} \times \frac{[\text{Standard}]}{2} \times n \quad (\mu\text{M})$$

where  $R_{\text{SAMPLE}}$ ,  $R_{\text{BLANK}}$  and  $R_{\text{STANDARD}}$  are OD or fluorescence readings of the Sample, Sample Blank and the Sample plus Standard respectively.  $n$  is the sample dilution factor. Note: The volume of the internal standard is 2× lower than the sample volume; thus, the internal standard concentration should be divided by 2. If the calculated maltose concentration is >500  $\mu\text{M}$  for the colorimetric assay, or >50  $\mu\text{M}$  for the fluorometric assay, dilute sample in dH<sub>2</sub>O and repeat assay. Multiply result by the dilution factor  $n$ .

Conversions: 1 mM maltose equals 34.23 mg/dL, or 342.3 ppm.

## Sample Data



Version: V.08.09.2018

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)