



Glucose Uptake Assay Kit (Fluorometric)

LS-K336-100 (100 Tests) • Store at -20°C

Introduction

Glucose uptake has a variety of methods and transporters, and depends upon the metabolic demand of the cell type and availability of glucose. There are over ten different facilitated diffusion glucose transporters which transport glucose down its concentration gradient without ATP hydrolysis. In the kidneys, secondary active transport is used to uptake glucose against its concentration gradient to ensure that very little glucose is excreted in urine. LSBio's fluorescent cell-based glucose uptake assay uses 2-deoxyglucose (2-DG), a widely used glucose analog because it can be taken up by glucose transporters and metabolized by endogenous hexokinase into 2-deoxyglucose 6-phosphate (2-DG6P). 2-DG6P accumulates intracellularly because it is not a suitable substrate for phosphoglucose isomerase, the next step in glycolysis. The cells are lysed, and excess NADP and glucose 6-phosphate dehydrogenase (G6PDH) is added to metabolize 2-DG6P and generate a molar equivalent amount of NADPH. The NADPH is then measured using a G6PDH recycling reaction to amplify the signal and generate a fluorescent signal measurable at $\lambda_{ex/em} = 530/585$ nm proportional to the concentration of 2-DG6P.



Key Features

- Safe. No radioactive material is used.
- Sensitive and Accurate. Detection limit of 0.1 μM and linearity up to 5 μM 2-DG6P.
- Simple and Convenient. Can be automated as a medium throughput assay for glucose transport in cells.

Applications

- Determination of glucose uptake in whole cells.
- Evaluation of effects of ligands or drugs on glucose transport.

FOR RESEARCH USE ONLY! Not for use in humans.

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

www.LSBio.com • (206) 464-1554 • TechnicalSupport@LSBio.com

Glucose Uptake Assay Kit (Fluorometric)

LS-K336-100 (100 Tests) • Store at -20°C



Components

Component	K336-100
	100 Tests
Assay Buffer	10 mL
G6P Reagent	1.5 mL
Enzyme A	120 µL
Enzyme B	120 µL
G6PDH Enzyme	120 µL
NADP	120 µL
2-DG Substrate	1.2 mL
Probe	750 µL
2-DG6P Standard	120 µL
NADP/NADPH Extraction Buffers	12 mL each

Materials Not Supplied

Triton X-100 (Sigma, cat # T8787); Phosphate Buffered Saline (Sigma, cat# P4417; can also be made yourself if desired); black cell culture 96-well plate: available separately at Sigma (CLS3603); deionized or distilled water; pipetting devices; cell culture incubators; centrifuge tubes; fluorescence plate reader capable of reading at $\lambda_{ex/em} = 530/585$ nm.

Storage

Store all reagents at -20°C. Shelf life of 6 months after receipt.

Assay Procedure

Important

1. To avoid cross-contamination, change pipette tips between additions of each reagent or sample. We recommend the use of a multi-channel pipette. Use separate reservoirs for each reagent.
2. It is recommended that samples be assayed in triplicate or higher.

Procedure

A. Culture, Starve and Treat Cells

1. Seed 100 µL of $1-10 \times 10^3$ adherent cells (or $1-5 \times 10^4$ suspension cells) into each well of a 96-well culture plate. Incubate for 4 hours or overnight at 37°C in a cell culture incubator.

Note: The cell number to be used depends on cell size and metabolic demand of glucose.

2. Incubate the cells with serum-less media for 4 hours or overnight to increase their glucose demand.
3. Starve the cells in glucose-less and serum-less media for 40 minutes. Add any drugs or experimental treatments to the starvation media at this step if desired. Make sure to include a control group without any experimental conditions.

FOR RESEARCH USE ONLY! Not for use in humans.

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

www.LSBio.com • (206) 464-1554 • TechnicalSupport@LSBio.com

Glucose Uptake Assay Kit (Fluorometric)

LS-K336-100 (100 Tests) • Store at -20°C



B. Add 2-Deoxyglucose

1. Add 10 μL of 2-DG substrate to each well. Incubate for 20 minutes or desired time.
2. Remove the media. Then wash the cells 3 times with 150 μL of ice cold PBS to remove excess 2-DG. Each wash should be performed for 30 seconds without shaking, try not to disturb the cells.

C. Lyse and Extract 2-Deoxyglucose 6-phosphate

1. Prepare Lysis buffer by adding 1% Triton X-100 to NADP Extraction Buffer. Prepare enough for 55 μL per well.
2. Standard Curve. Prepare 5 μM 2-DG6P Premix by mixing 5 μL 5 mM Standard and 4995 μL distilled water. Dilute standard as follows.

No	Premix + H ₂ O	2-DG6P (μM)
1	100 μL + 0 μL	5.0
2	60 μL + 40 μL	3.0
3	30 μL + 70 μL	1.5
4	0 μL + 100 μL	0

Transfer 50 μL standards into separate wells of the cell plate.

3. Transfer 50 μL of lysis buffer to each well with cells, place the plate on a rotary shaker for 5 minutes, and then incubate the plate at 80°C for 10 minutes.
4. Add 50 μL of NADPH extraction buffer to each well with cells, and 50 μL of dH₂O to the standard wells. Cool the plate in a -20°C freezer for 5 minutes followed by 10 min on the bench top (alternatively you may also cool the plate to room temperature on the bench top (~30 min)).
5. Working Reagent 1 (WR1) Preparation. For each reaction well, prepare WR1 by mixing 10 μL Assay Buffer, 1 μL G6PDH Enzyme, 1 μL NADP. Add 10 μL of WR1 to all wells and incubate at 37°C for 60 minutes.

D. Extract and Measure NADPH

1. Add 50 μL of NADPH extraction buffer to all wells and incubate at 80°C for 15 minutes.
2. Add 50 μL of NADP extraction buffer to all wells, and cool the plate in a -20°C freezer for 5 minutes (alternatively you may also leave the plate in a refrigerator or on the bench top).
3. Transfer 50 μL of sample and standard from each well in the cell plate into separate wells in a black 96-well plate.
4. Working Reagent Preparation. For each reaction well, prepare Working Reagent by mixing 45 μL Assay Buffer, 1 μL Enzyme A, 1 μL Enzyme B, 10 μL G6P Reagent and 5 μL Probe. Fresh reconstitution is recommended.
5. Transfer 50 μL of Working Reagent into each well. Read the plate at $\lambda_{\text{ex/em}} = 530/585 \text{ nm}$ for 20 minutes. Use data from time zero and 20 minutes (F_0) and (F_{20}).

Note: If fluorescent signal for any sample is higher than the fluorescence of the 5 μM standard, dilute the sample in dH₂O and repeat Steps 3-5 in D. Extract and Measure NADPH. Multiply the results by the dilution factor.

FOR RESEARCH USE ONLY! Not for use in humans.

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

www.LSBio.com • (206) 464-1554 • TechnicalSupport@LSBio.com

Glucose Uptake Assay Kit (Fluorometric)

LS-K336-100 (100 Tests) • Store at -20°C


Calculations

First compute the ΔF for each standard and sample by subtracting F_0 from F_{20} . If duplicate or triplicate samples were performed, calculate the mean ΔF intensities for the Sample wells. Plot the standard ΔF 's and determine the slope. The concentration of 2-DG6P is calculated as follows:

$$[2\text{-DG6P}] = \frac{\Delta F_{\text{Sample}} - \Delta F_{\text{Blank}}}{\text{Slope}} \times n \quad (\mu\text{M})$$

where ΔF_{Sample} is the mean ΔF of the samples, and ΔF_{Blank} is the ΔF of the water, standard #4 blank. Slope is the slope of the standard curve and n is the dilution factor.

Sample Data



Glucose Transport Assay in PANC-1 Cells. Cells were seeded, starved, and treated according to protocol. Left: PANC-1 cell titration in the absence and presence of 1 mM phloretin. Right: glucose transport inhibition curve with phloretin. PANC-1 cells were seeded at 10,000 cells per well.

Version: V.08.09.2018

FOR RESEARCH USE ONLY! Not for use in humans.

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

www.LSBio.com • (206) 464-1554 • TechnicalSupport@LSBio.com