


# Pyruvate Dehydrogenase (PDH) Activity Assay Kit (Colorimetric)

LS-K619-100 (100 Tests) • See Storage Conditions Below



## Introduction

Pyruvate Dehydrogenase (PDH) (EC 1.2.4.1) has a vital role in carbohydrate metabolism. It forms a well-characterized enzyme complex with dihydrolipoyl transacetylase (E2) and dihydrolipoyl dehydrogenase (E3). PDH converts pyruvate into acetyl-CoA in the presence of NAD and CoA, and links glycolysis to the citric acid cycle. PDH activity is inhibited by high intracellular ratios of ATP/ADP, NADH/NAD or Acetyl-CoA/CoA. In humans, PDH deficiency reduces mitochondrial function and is linked to neurodegenerative diseases. PDH deficiency is X-linked; it results in 2 forms of abnormality: a metabolic form (lactic acidosis) and a neurological form (seizure and/or neuropathological spasm). Recent studies show that PDH is a target of oncogene-induced senescence; activation of PDH enhances pyruvate utilization and increases respiration and redox stress. LSBio's PDH Assay Kit provides a quick and easy way for monitoring PDH activity in various samples. In the assay, PDH converts pyruvate into an intermediate, which reduces the developer to a colored product with strong absorbance at 450 nm. The assay is simple, sensitive and can detect pyruvate dehydrogenase activity lower than 0.1 mU in a variety of samples.



## Applications

- Measurement of pyruvate dehydrogenase activity in various tissues/cells
- Analysis of cell signaling pathway

## Sample Types

- Animal tissues: heart, liver, muscle, etc.
- Purified mitochondria
- Cell culture: adherent or suspension cells

## Components

Component	K619-100	Cap Code
	100 Tests	
PDH Assay Buffer	25 ml	WM
PDH Substrate (Lyophilized)	1 vial	Blue
PDH Developer (Lyophilized)	1 vial	Red
NADH Standard (Lyophilized)	1 vial	Yellow
PDH Positive Control	10 µl	Orange

## Materials Not Supplied

- 96-well clear plate with flat bottom
- Mitochondria Isolation Kit (LSBio, Cat. # LS-K158)
- Saturated ammonium sulfate

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)

# Pyruvate Dehydrogenase (PDH) Activity Assay Kit (Colorimetric)

LS-K619-100 (100 Tests) • See Storage Conditions Below



---

## Storage Conditions and Reagents Preparation

Store kit at -20°C, protected from light. Briefly centrifuge small vials prior to opening. Read the entire protocol before performing the assay.

- PDH Assay Buffer: Warm to room temperature before use. Store at either 4°C or -20°C.
- PDH Substrate: Reconstitute with 220 µl dH<sub>2</sub>O. Store at -20°C. Keep on ice while in use. Use within two months.
- PDH Developer: Reconstitute with 220 µl dH<sub>2</sub>O. Gently pipette up and down to dissolve completely. Store at -20°C. Use within two months.
- NADH Standard: Reconstitute with 400 µl dH<sub>2</sub>O to generate 1.25 mM NADH Standard solution. Aliquot and store at -20°C. Keep on ice while in use. Use within two months.
- PDH Positive Control: Add 100 µl PDH Assay Buffer to the Positive Control and mix thoroughly. Aliquot and store at -20°C. Keep on ice while in use. Use within two months.

## Assay Procedure

1. Sample Preparation: Rapidly homogenize tissue (10 mg) or cells ( $1 \times 10^6$ ) with 100 µl ice-cold PDH Assay Buffer, and keep on ice for 10 min. Centrifuge at 10,000 x g for 5 min. and transfer the supernatant to a fresh tube. Add 5-50 µl sample per well & adjust the volume to 50 µl with PDH Assay Buffer. To check PDH activity in mitochondria, isolate the mitochondria from fresh tissue or cells using LSBio's Mitochondria Isolation Kit for Tissue and Cultured Cells (LSBio, Cat. # LS-K158). Add 5-50 µl of isolated mitochondria per well, adjust the volume to 50 µl with PDH Assay Buffer. For the PDH positive control, take 1-10 µl of PDH Positive Control into desired well(s) and adjust the final volume to 50 µl with PDH Assay Buffer.

### Notes:

- a) For unknown samples, we suggest testing several doses to ensure the readings are within the Standard Curve range.
  - b) For samples exhibiting significant background, prepare parallel sample well(s) as background controls.
  - c) Small molecules in some tissues such as liver may interfere with the assay. To remove small molecules, we suggest using an ammonium sulfate method. Pipette 50-100 µl of lysate into a fresh tube, add 2X volume of saturated ammonium sulfate (4.1 M, room temperature) and keep on ice for 20 min. Spin down at 10,000 X g for 5 min., carefully remove and discard the supernatant, and re-suspend the pellet to the original volume with PDH Assay Buffer.
2. NADH Standard Curve: Add 0, 2, 4, 6, 8 and 10 µl of 1.25 mM NADH Standard into a series of wells in 96-well plate to generate 0, 2.5, 5.0, 7.5, 10 and 12.5 nmol/well of NADH Standard. Adjust the volume to 50 µl/well with PDH Assay Buffer.

---

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121  
[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)

# Pyruvate Dehydrogenase (PDH) Activity Assay Kit (Colorimetric)

LS-K619-100 (100 Tests) • See Storage Conditions Below

- Reaction Mix: Mix enough reagents for the number of assays to be performed. For each well, prepare 50  $\mu$ l Mix containing:

	Reaction Mix	*Background Control Mix
PDH Assay Buffer	46 $\mu$ l	48 $\mu$ l
PDH Developer	2 $\mu$ l	2 $\mu$ l
PDH Substrate	2 $\mu$ l	----

Mix and add 50  $\mu$ l of the Reaction Mix to each well containing the Standard, Positive Control and test samples.

\*For background correction, add 50  $\mu$ l of Background Control Mix (without substrate) to sample background control well(s) and mix well.

- Measurement: Measure absorbance immediately at 450 nm in kinetic mode for 10-60 min. at 37°C.

Note: Incubation time depends on the pyruvate dehydrogenase activity in samples. We recommend measuring the OD in kinetic mode, and choosing two time points ( $T_1$  &  $T_2$ ) in the linear range to calculate the pyruvate dehydrogenase activity of the samples. The NADH Standard Curve can be read in Endpoint mode (i.e., at the end of the incubation time).

- Calculation: Subtract 0 Standard reading from all readings. Plot the NADH Standard Curve. If sample background control reading is significant, subtract the background control reading from its paired sample reading. Calculate the pyruvate dehydrogenase activity of the test sample:  $\Delta OD = A_2 - A_1$ . Apply the  $\Delta OD$  to the NADH Standard Curve to get B nmol of NADH generated during the reaction time ( $\Delta T = T_2 - T_1$ ).

$$\text{Sample Pyruvate Dehydrogenase Activity} = B / (\Delta T \times V) \times D = \text{nmol/min/ml} = \text{mU/ml}$$

Where: B = NADH amount from Standard Curve (nmol),  $\Delta T$  = reaction time (min.), V = sample volume added into the reaction well (ml), D = Dilution Factor.

Unit Definition: One unit of pyruvate dehydrogenase is the amount of enzyme that generates 1.0  $\mu$ mol of NADH per min. at pH 7.5 at 37°C.

## Sample Data




Figure: (a) NADH standard curve; (b) Pyruvate Dehydrogenase activity in Jurkat mitochondria (100  $\mu$ g); (c) Pyruvate Dehydrogenase specific activity was calculated in mitochondria prepared from mouse heart (10  $\mu$ g), *S. Cerevisiae* (10  $\mu$ g), MCF-7 cells (30  $\mu$ g) or Jurkat cells (35  $\mu$ g). Assays were performed following the kit protocol.

Version: V.08.09.2018

**FOR RESEARCH USE ONLY! Not for use in humans.**

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121

[www.LSBio.com](http://www.LSBio.com) • (206) 464-1554 • [TechnicalSupport@LSBio.com](mailto:TechnicalSupport@LSBio.com)