

α -Amylase Assay Kit (Colorimetric)

LS-K167-100 (100 Tests) • Store at -20°C



Introduction

AMYLASE belongs to the family of glycoside hydrolase enzymes that break down starch into glucose molecules by acting on α -1,4-glycosidic bonds. The α -amylases (EC 3.2.1.1) cleave at random locations on the starch chain, ultimately yielding maltotriose and maltose, glucose and "limit dextrin" from amylose and amylopectin. In mammals, α -amylase is a major digestive enzyme. Increased enzyme levels in humans are associated with salivary trauma, mumps due to inflammation of the salivary glands, pancreatitis and renal failure.

Simple, direct and automation-ready procedures for measuring amylase activity are very desirable. This α -amylase assay Method involves two steps: (1). α -amylase in the sample hydrolyzes starch and the product is rapidly converted to glucose by α -glucosidase and hydrogen peroxide by glucose oxidase; (2). hydrogen peroxide concentration is determined with a colorimetric reagent.

Key Features

- Sensitive and accurate. Linear detection range 0.3 to 50 U/L α -amylase in 96-well plate assay.
- Convenient. The procedure involves adding a single working reagent, incubation for 15 min, followed by the detection reagent and a 20-min incubation and reading the optical density at 585 nm.

Applications

- Determination of α -amylase activity in blood, saliva, urine, grains and other agricultural samples.

Components

Component	K167-100
	100 Tests
Assay Buffer (pH 7.0)	20 mL
Substrate	120 μ L
Detection Reagent	20 mL
Enzyme A	120 μ L
Enzyme B	120 μ L
Glucose Standard	1 mL

Materials Not Supplied

Pipetting devices, centrifuge tubes, clear flat-bottom 96-well plates, plate reader, and optionally membrane filters (e.g. Microcon YM-10 from Millipore).

Storage

Kit is shipped on ice. Store all components at -20°C. Shelf life: 6 months after receipt.

FOR RESEARCH USE ONLY! Not for use in humans.

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121
www.LSBio.com • (206) 464-1554 • TechnicalSupport@LSBio.com

α -Amylase Assay Kit (Colorimetric)

LS-K167-100 (100 Tests) • Store at -20°C



Assay Procedure

Reagents. Equilibrate all components to room temperature. Keep thawed Enzyme Mix in a refrigerator or on ice. The substrate may have precipitates. Prior to use, vortex tube to dissolve precipitates; gentle swirl the Detection Reagent bottle.

Sample preparation. Ideally samples are assayed fresh. When stored frozen, α -amylase is stable for one month. Ascorbic acid, heparin, EDTA, EGTA, citrate, SDS, Tris (> 8mM) and ethanol (>0.4%) interfere and should be avoided in sample preparation. If glucose is present in the sample, treat the samples as described in General Considerations. It is prudent to perform a pilot test with samples at various dilutions. Recommended dilution: serum 50-fold, saliva 2,000-fold in Assay Buffer prior to assay.

1. Prepare 400 μ M Glucose Standard by mixing 10 μ L of the provided (300 mg/dL) standard with 406 μ L Assay Buffer. Transfer 10 μ L Assay Buffer, 10 μ L 400 μ M glucose, and 10 μ L of each sample into separate wells of a clear flat-bottom 96-well plate.
2. Prepare enough Working Reagent for each well by mixing 40 μ L Assay Buffer, 0.5 μ L Substrate, 1 μ L Enzyme A, 1 μ L Enzyme B.
Transfer 40 μ L Working Reagent to each well. Incubate for 15 min at room temperature (25°C).
3. Add 150 μ L Detection Reagent to each well. Mix and incubate for 20 min at room temperature (25°C). Read OD585nm (540-610nm) on a plate reader.

Calculations

The Amylase activity is calculated as

$$\text{Activity} = \frac{\text{OD}_{\text{SAMPLE}} - \text{OD}_{\text{BUFFER}}}{\text{OD}_{\text{STD}} - \text{OD}_{\text{BUFFER}}} \times \frac{400}{t \text{ (min)}} \times n \quad (\text{U/L})$$

$\text{OD}_{\text{SAMPLE}}$, OD_{STD} and $\text{OD}_{\text{BUFFER}}$ are optical density values of the sample, the 400 μ M glucose standard and Assay Buffer. t is the incubation time. $t = 15$ min in the standard protocol. n is the dilution factor ($n = 50$ for serum, 2000 for saliva). One unit of enzyme catalyzes the production of 1 μ mole of glucose per min under the assay conditions.

Note: if the calculated activity is higher than 50 U/L, dilute sample in Assay Buffer and repeat assay. Multiply the results by the dilution factor.

General Considerations

For samples known to contain glucose, use a membrane filter (e.g. Microcon YM-10 from Millipore) to remove glucose: load 50 μ L sample in a Microcon YM-10 (10 kDa cutoff) and add 500 μ L Assay Buffer. Centrifuge at 14000 rpm for 30 min, check level of sample, ideally the sample level will be less than 50 μ L. Add 500 μ L Assay Buffer and repeat the centrifugation. Measure final sample volume with a pipetman and calculate dilution factor $n = \text{final sample volume}/50 \mu\text{L}$.


FOR RESEARCH USE ONLY! Not for use in humans.

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121
www.LSBio.com • (206) 464-1554 • TechnicalSupport@LSBio.com

α -Amylase Assay Kit (Colorimetric)

LS-K167-100 (100 Tests) • Store at -20°C

Sample Data



Standard Curve in 96-well plate assay

Version: V.08.09.2018

FOR RESEARCH USE ONLY! Not for use in humans.

LifeSpan BioSciences, Inc. • 2401 Fourth Avenue, Suite 900, Seattle, WA 98121
www.LSBio.com • (206) 464-1554 • TechnicalSupport@LSBio.com